

# Colour Eye Tests


## Ishihara Colour Chart

### Testing colour perception

Ishihara colour charts are ideal for the screening of hereditary red-green colour vision defects. Available in four different versions, book format.

Testing distance: 75 cm

- 47640 Ishihara's tests for colour deficiency, 38 plates edition
- 47600 Ishihara's tests for colour deficiency, 24 plates edition
- 47610 Ishihara's tests for colour deficiency, 14 plates edition
- 47620 Ishihara's design charts for colour deficiency of unlettered persons, 10 plates




## Matsubara Colour Charts

### Colour vision test plates for the infants

Matsubara colour charts are ideal for the screening of hereditary red-green colour vision defects. 10 test plates in book format, suitable from pre-school age up.

Testing distance: 75 cm

- 47625 Matsubara colour test plates for the infants, book format, 10 plates, incl. 6 matching charts



## Colour Arrangement Tests

Quick screening test for differentiation between mild and pronounced colour vision defects.

- 47650 Farnsworth colour test panel 15 D (saturated), in wooden box, incl. evaluation form
- 47651 Lanthony colour 15 HUE (desaturated), in wooden box, incl. evaluation form

